

Fifteen years of real-time stability data at room temperature and accelerated aging: Validation of DNAsell[®] encapsulation for sustainable DNA data storage

Marthe COLOTTE, Aurélie LUIS, Delphine COUDY, Magali MILHAU, Sophie TUFFET, Jacques BONNET

Marthe COLOTTE, PhD
Imagine CTO
colotte@imagine.eu

DNA data storage

- Will contribute to resolve the data volume explosion crisis
- Parameters to take into account / challenges to address
 - Scale of storage
 - Information density
 - DNA stability* / data retrieval




Stability often overlooked

- DNA degradation mechanisms and factors
- DNAsell® / RNAsell® technology
- Evaluation of stability
 - Accelerated aging (Arrhenius model)
 - Real time aging @RT


*We'll talk about data encoded in DNA sequence, not in its structure

Main mechanisms of degradation



Consequences


- The main and most deleterious event:
chain breakage
- The relative importance of hydrolysis and oxidation vary according to:
 - **Water content**
 - Environment (contaminants, co-extractants, ...)



Bonnet *et al.* Nucleic acids research 2010
10.1093/nar/gkp1060 (corrected)

Water is the enemy

How hermetic is a container/sealing?




- Counter intuitive results!
- Inner RH will gradually equilibrate with outer RH
- In a closed but non hermetic system, water uptake has to be taken into account during stability testing (modeling can be tricky)

Bonnet *et al.* Nucleic acids research 2010
10.1093/nar/gkp1060

New container highly needed!


The Imagen technology

- Imagen's technology: completely remove water, oxygen and ozone and maintain these conditions for efficient and long-term storage.
- Imagen stainless steel capsules: DNAShell® and RNAShell®
 - Stainless steel shell
 - Glass vial
 - Stainless steel cap
- The capsule is sealed by YAG laser welding allowing to maintain the desiccated sample under an inert atmosphere.
- A unique 2D barcodes is engraved on each capsule for full traceability.
- Target markets: biobanking, molecular diagnostics, DNA data storage





Monitoring chain breakage with pDNA

Only one chain break on one of its two strands is sufficient to relax a supercoiled plasmid. The supercoiled plasmid content (SCC) is thus the fraction of intact DNA in the sample. The relaxation is a first order reaction (exponential).



Monitoring chain breakage with gDNA



L_{MAX} = size of the maximum of intensity on the gel (single strand)

- L_{MAX} : from k_T and time t (estimation after Arrhenius equation)
- k_T from time t and L_{MAX} (from gel)

Long term stability evaluation

■ Accelerated aging

- Conduct degradation kinetics at several high temperatures T
- Measure the degradation rates k_T (chain breaks/nt /s)
- Plot $-\log k_T$ vs $1/T$
 - => Arrhenius equation
 - => extrapolation of degradation rates @RT (or any temperature)

*Bonnet, J., *et al.*, Nucleic acids research 2010 10.1093/nar/gkp1060

*Coudy, D., *et al.*, Long term conservation of DNA at ambient temperature. Implications for DNA data storage. PLoS One, 2021. 16(11)

* Fabre, A-L. *et al.*, An efficient method for long-term room temperature storage of RNA. Eur J Hum Genet 2014 Mar;22(3):379-85

* Fabre, A-L. *et al.*, High DNA stability in white blood cells and buffy coat lysates stored at ambient temperature under anoxic and anhydrous atmosphere PLoSOne 2017 10.1371/journal.pone.0188547

Long term stability evaluation

- Real time evaluation of degradation rate at room temperature
 - To verify the accelerated aging results
 - To take potential unforeseen parameters into account

* Colotte *et al.* Biopreserv Biobank. 2011 10.1089/bio.2010.0028

This work (unpublished)

* Fabre, A-L. *et al.*, High DNA stability in white blood cells and buffy coat lysates stored at ambient temperature under anoxic and anhydrous atmosphere PLOSOne 2017 10.1371/journal.pone.0188547

Material and methods

■ Source biospecimens and extraction methods

Biological material	Purified		In biospecimen	
	DNA	RNA	DNA	RNA
Cell lines	Organic Salt precip Silica	Phenol/Chlo Trizol	-	Fixed cells Cells on paper
Bacteria	Organic Salt precip Silica Ion exchange	Phenol/Chlo Silica Trizol	-	-
Rat liver	Organic Salt precip	Phenol/Chlo Silica Trizol	-	-
Blood	Organic Salt precip	Paxgene	Lysed WBC Lysed BC	-
Saliva	Organic Salt precip	-	-	-
Plant	Organic Silica	-	-	-

■ Processing

- Collection of the biospecimens
- Purified NA: extraction + QC
- NA in biospecimen: lysis for DNA in blood cells / Cell fixation or deposition on paper for RNA
- Addition of stabilization proprietary solutions
- Aliquoting in the minicapsules
- Desiccation
- Encapsulation (laser sealing under argon)

Material and methods

- Accelerated aging – simulated 25-100 yrs: no degradation (not shown)

- Real-time long term storage



← Capsules stored at room temperature without moisture and temperature control (15°C - 30°C / 50 % RH)



← Control capsules stored at -20°C

← Controls in solution stored at -20°C (DNA) or -80°C (RNA)




- Post-storage analysis

- Opening of the minicapules with a shellOpener
- Rehydration
- Extraction for biospecimen (salt precipitation for DNA / Trizol for RNA)
- QC (not shown)
- Electrophoresis agarose gels/ Bioanalyzer (RIN)
- RT-qPCR / RT-dPCR analysis (not shown)

Results

Purified DNA






Heat denaturation + 0.8% agarose gel – TAE buffer – EtBr staining – Size in kb

$$L_{\max} > 20 \text{ kb after } t = 14.5 \text{ yrs} \Rightarrow k_{25^\circ\text{C}} < 1.09 \times 10^{-13} \text{ /nt/s}$$

No difference between the frozen controls and the capsules (same as initial sizes, not shown).


Plasmid DNA



$$m/m_0 = 0.93 \text{ after } t = 15 \text{ yrs} \Rightarrow k_{25^\circ\text{C}} = 1.4 \times 10^{-14} \text{ /nt/s}$$



Colotte *et al.* Biopreserv Biobank. 2011
10.1089/bio.2010.0028

Accelerated vs real-time




Redrawn from Coudy, D., *et al.*, Long term conservation of DNA at ambient temperature. Implications for DNA data storage. PLoS One, 2021. 16(11)


Purified RNA



DNA and RNA in biospecimen



*Heat denaturation + 0.8% agarose gel –
TAE buffer – EtBr staining - Size in kb*



Conclusions

- DNAsell /RNAsell is the most advanced available technology for DNA/RNA stabilization, with 15 years of recorded performance
- To build confidence in DNA data storage, it is imperative to:
 - Use validated methods (SNIA standard) & run all necessary controls
 - Use high-sensitive methods to demonstrate stability (easy to hide poor performance)



imagine

Smart Biopreservation

Thank you !

Marthe COLOTTE, PhD
Imagine CTO
colotte@imagine.eu

Publications for DNA preservation

External

In house

■ DNA



Stabilité chimique et conformationnelle de l'ADN à l'état sec et à température ambiante

Thèse de Marthe Colotte (2008)

Chain and conformation stability of solid-state DNA: implications for room temperature storage

DOI: [10.1093/nar/gkp1060](https://doi.org/10.1093/nar/gkp1060)

Simultaneous assessment of average fragment size and amount in minute samples of degraded DNA - DOI: [10.1016/j.ab.2009.02.003](https://doi.org/10.1016/j.ab.2009.02.003)

Adverse effect of air exposure on the stability of DNA stored at room temperature

DOI: [10.1089/bio.2010.0028](https://doi.org/10.1089/bio.2010.0028)

Novel procedure for high yield recovery of traces amounts of DNA stored at room temperature

Poster - ESBB

Assessment of DNA encapsulation, a new room-temperature DNA storage method

DOI: [10.1089/bio.2013.0082](https://doi.org/10.1089/bio.2013.0082)



Evaluation of DNA/RNAs shells for room temperature nucleic acids storage

DOI: [10.1089/bio.2014.0060](https://doi.org/10.1089/bio.2014.0060)



Preservation of biospecimens at ambient temperature: special focus on nucleic acids and opportunities for the biobanking community

DOI: [10.1089/bio.2015.0022](https://doi.org/10.1089/bio.2015.0022)

Quality Matters: 2016 Annual Conference of the National Infrastructures for Biobanking

DOI: [10.1089/bio.2016.0053](https://doi.org/10.1089/bio.2016.0053)

Ensuring the Safety and Security of Frozen Lung Cancer Tissue Collections through the Encapsulation of Dried DNA

DOI: [10.3390/cancers10060195](https://doi.org/10.3390/cancers10060195)



DNAsell Protects DNA Stored at Room Temperature for Downstream Next-Generation Sequencing Studies

DOI: [10.1089/bio.2018.0129](https://doi.org/10.1089/bio.2018.0129)



■ DNA reference material

A novel method for room temperature distribution and conservation of RNA and DNA reference materials for guaranteeing performance of molecular diagnostics in onco-hematology: A GBMHM study

DOI: [10.1016/j.clinbiochem.2015.04.004](https://doi.org/10.1016/j.clinbiochem.2015.04.004)



■ White Blood Cells & Buffy coat

High DNA stability in white blood cells and buffy coat lysates stored at ambient temperature under anoxic and anhydrous atmosphere

DOI: [10.1371/journal.pone.0188547](https://doi.org/10.1371/journal.pone.0188547)

■ Synthetic DNA

→ DNA data storage

Long term conservation of DNA at ambient temperature. Implications for DNA data storage

DOI: [10.1371/journal.pone.0259868](https://doi.org/10.1371/journal.pone.0259868)

An Empirical Comparison of Preservation Methods for Synthetic DNA Data Storage

DOI: [10.1002/smt.202001094](https://doi.org/10.1002/smt.202001094)

UNIVERSITY of WASHINGTON



Microsoft Research

Publications for RNA preservation

■ RNA

2014

An efficient method for long-term room temperature storage of RNA

DOI: [10.1038/ejhg.2013.145](https://doi.org/10.1038/ejhg.2013.145)

2014

Evaluation of DNA/RNAs shells for room temperature nucleic acids storage

DOI: [10.1089/bio.2014.0060](https://doi.org/10.1089/bio.2014.0060)

华大基因
BGI

2017

Long-term room temperature storage of dry ribonucleic acid for use in RNA-Seq analysis

DOI: [10.1089/bio.2017.0024](https://doi.org/10.1089/bio.2017.0024)

华大基因
BGI

■ RNA Reference material

2015

A novel method for room temperature distribution and conservation of RNA and DNA reference materials for guaranteeing performance of molecular diagnostics in onco-hematology: a GBMHM study

DOI: [10.1016/j.clinbiochem.2015.04.004](https://doi.org/10.1016/j.clinbiochem.2015.04.004)



■ Synthetic RNA

→ Reference material SARS-CoV-2

2023

Reference materials for SARS-CoV-2 molecular diagnostic quality control: validation of encapsulated synthetic RNAs for room temperature storage and shipping

DOI: [10.1101/2023.08.28.555008](https://doi.org/10.1101/2023.08.28.555008)



→ Certified Reference material SARS-CoV-2

2024

Certification of the identity and the copy number concentration of synthetic single-stranded RNA including fragments of the SARS-CoV-2 genome and part of the human RNase P gene: EURM[®]-014k – JRC Reference Material report.



Publications for Biospecimen

■ Blood biomarkers (for diagnostics)

2017

Stability of newborn screening markers in dried-blood spot (DBS): the innovative imagene solution.

Poster - SFEIM



revvity



■ Live bacteria & viruses

2017

ANVBIS3 – Acides nucléiques, virus & bactéries d'intérêt en stockage standardisé et sécurisé

Poster - 6ème Forum DGA innovation



■ Spermatozoa (*for nucleus transfert*)

2023

Reviving Vacuum-Dried Encapsulated Ram Spermatozoa via ICSI after 2 Years of Storage

DOI: manuscript accepted



2025

Fertility preservation of vacuum-dried ram spermatozoa stored for four years at room temperature

Theriogenology – Vol 239, June 2025, 117390

